

# XPERIMATE QUADRA TOWER

## LIQUID HANDLING WORKSTATION



SYNCHRON

Laboratory automation & engineering

- Easy to use
- 96 Individual pipettors can be 24 well
- 450 $\mu$ L to less than 10 $\mu$ L pipetting range
- 2% CV at 10 $\mu$ L
- Manual push-button operation
- Fully Programmable
- No external computer required
- One year warranty parts & labor

## MULTIWELL PIPETTING



### DESCRIPTION OF OPERATION

The XPERIMATE QUADRA TOWER is designed to fill the gap between 8 and 12-channel, hand-held pipettors and high-level automation. It may be used in a simple, manual mode. Individual push buttons provide the basic functions of aspirating and dispensing by simply entering the volume desired. Other push buttons select tip-loading and tip-shucking. Where feasible, the tips may be washed for economy of operation. More sophisticated pipetting sequences can be fully programmed in a simple, intuitive manner and selected for repetitive applications. This is accomplished from the front control panel. It can also be integrated into a robotic system. In essence, the mode of operation is designed for operational convenience by a variety of users, using different skill levels.

Two pipetting stations accommodate liquid transfers from reservoirs or deepwell microplates. Two other stations accommodate 96 well microplates. The Quadra tip design is unique. The maximum volume is 450 $\mu$ L, particularly useful if aliquoting to repetitive plates. The long, narrow end will access the bottom of a deepwell plate. The design also permits the pipeline pipetting of individual reagents using air-gap separation. The unit has a CV of 2% at 10 $\mu$ L.

### AVAILABLE ACCESSORIES

A peristaltic pump is available to continuously change the water in the tip washing station. An ultrasonic tip washing station provides additional economy of operation for more demanding applications. The 96 well pipettor head assembly is standard for the Tower. However, it can also be supplied with small volume, positive displacement pipettor heads in either 96 well or 384 well. These pipettors can provide precision of 2-3% at the 0.5 $\mu$ L level.



For more information:

[WWW.SYNCHRONLAB.COM](http://WWW.SYNCHRONLAB.COM)